

Abdomen/Pelvis 6.21
Low Dose Follow Up Renal stone
Fixed MA
VCT/LS 16

Indications: Follow up of known urinary calculus

Patient Preparation: Full bladder. Clamp Foley.
Oral contrast: None.
IV: Contrast: None.

****Position patient PRONE****

Group 1: Non-contrast – Top of kidneys to bottom of symphysis pubis

Scan Mode	Thickness Speed Pitch	Table Interval	SFOV	kV	MA	Prep Time (sec)
0.8 sec Helical Full	5 27.50 1.375:1	5	Large	120	75	None

Algorithm: Standard

Reorient image to standard supine position.

Network: Group 1 (5mm) to PACS. Group 2 (1.25mm) auto transmits AWBR.

Notes: This is a fixed MA (60 MAS)

Written: 01/10/2010