

Routine Head Trauma 1.2

LightSpeed16

Indication: Trauma

Patient Preparation: None

Group 1:

Scan Mode	Thickness Speed Pitch	Table Interval	SFOV	kV	Auto mA/ Noise Index	Prep Time (sec)	Recon Type
1.0 sec Axial Full	5 2i	5	Head	120	380 8	N	Bone Standard

Algorithm: Standard and bone. Show recon 1 Standard 2 Bone 3 Standard

Network: Recon 1 and 2 to PACS.

07/06