

EXAM _____ PULMONARY _____

COIL _____ TORSO PA _____

PLANE	MODE	PULSE SEQ.	IMAGE OPTIONS	TE	TR	FLIP ANGLE	ECHOS	FOV	SLICE THICKNESS INTERSPACE	MATRIX	NEX	FREQ. DIR.
COR LOC	2D	SSFSE	FAST,SS	MIN	MIN	_	62.50	48	8-2	256X128	1	S/I
AX LOC	2D	SSFSE	FAST,SS	MIN	MIN	_	62.50	48	8-2	256X128	1	R/L
TOF	3D	Vasc TOF SPGR	VBW,FAST,ZIP2 ,ZIP512,Smart Prep	MIN	_	45	31.25	32	2.6/32	256X128	1	UNS WAP

Comments: USE CONTRAST INJECTOR

Written: 05/2007
 Revised: 01/2009