

Weight Limits  
Radiology Equipment

Dept	Equipment	Weight limit	Restrictions
DIAG	Siemens Sireskop's (ACC)	400	do not move table top on anyone over 300
DIAG	GE Flouro (2nd floor)	350	do not move table top on anyone over 300
MRI	GE Signa Horizons	350	
MRI	GE Signa Horizons	350	
CT	GE LightSpeed 16 (FWCC/CT)	400	with .25mm positional accuracy or 450lbs with 1mm positional accuracy
CT	GE VCT 64 (CT)	500	
CT	Siemens (ER)	660	
Rad Onc	Philips CT ACQSIM	450	85 cm bore
Rad Onc	Elekta Linear Accelerator	450	Max. weight at <b>end</b> of table 220 lbs
Nuclear	GE DST-XLI	400	
Nuclear	Siemens E.cam	300	
Nuclear	Siemens Symbia S	400	
Special	GE LCA-Table Limit	350	Do not Extend Table and Perform CPR with Max LBS
Special	Philips Allura Xper FD 20 10	450	Do not Move Table over 450LBS