

Thoracic Spine 5.17

Helical LightSpeed16

Indication: Trauma or Routine

Patient Preparation: None

Group 1: C7 to L1

Scan Mode	Thickness Speed Pitch	Table Interval	SFOV	kV	Auto mA/ Noise Index	Prep Time (sec)	Recon Type
1.0 sec Helical Full	2.5 5.62 0.562:1	5	Large 16cm	120	380 7	N	Bone Standard

Algorithm: Standard and bone. Show recon 1 Standard 2 Bone 3 25cm Bone

Reformatting: Sagittal curved coronal off 25cm bone

Network: Recon 1, 2 and reformatted images to PACS.