

EXAM__PEDS MRCP_____

COIL HEAD/QUADKNEE

PLANE	MODE	PULSE SEQ.	IMAGE OPTIONS	TE	TR	FLIP ANGLE	ECHOS	FOV	SLICE THICKNESS INTERSPACE	MATRIX	NEX	FREQ. DIR.
3 PL	2D	LOC	-----	---	---	-----	-----	24	5-5	256X128	1	-----
COR SSFSE	2D	SSFSE	VBW,Fast,SS	80	2000	31.25	-----	18	4-1	256X192	1	UN SWAP
COR T2	2D	FSE-XL	FC,NPW,VBW Fast, EDR	102	2400	15.63	12	18	4-1	256X224	1	UN SWAP
AX T2	2D	SSFSE	FC,VBW,EDR Fast, SS	80	2000	31.25	-----	18	3-0.3	256X224	1	UN SWAP
MRCP THICK	2D	SSFSE	VBW, Fast,SS	767	4000	31.25	-----	18	20.0-0	256X224	----	UN SWAP
MRCP THIN	2D	SSFSE	VBW, Fast,SS	80	4000	31.25	-----	18	3-0	256X224	-----	UN SWAP

Comments: _____