

EXAM__PRE ARTHROGRAM_____

COIL _SHOULDER PA_____

| PLANE | MODE | PULSE SEQ. | IMAGE OPTIONS | TE | TR | FLIP ANGLE | ECHOS | FOV | SLICE THICKNESS INTERSPACE | MATRIX | NEX | FREQ. DIR. |
|---------------|------|------------|-------------------------------|-----|------|------------|----------|-----|----------------------------|---------|-----|----------------|
| 3PL LOC | 3PL | LOC | NPW | - | - | - | - | 24 | 5-2 | 256X128 | 2 | - |
| COR T2 | 2D | FSE-XL | NPW,VBW,EDR, Fast,SCIC ZIP512 | 102 | 2750 | - | 16/15.63 | 16 | 4-0 | 256X192 | 2 | UNS WAP |
| COR T2 Fatsat | 2D | FSE-XL | Fast,FC,ZIP512 | 102 | 2200 | - | 4/20.83 | 16 | 4-0 | 256X192 | 2 | UNS WAP Fatsat |
| SAG T2 | 2D | FSE-XL | Fast,FC,ZIP512 | 85 | 2800 | - | 12/20.83 | 16 | 4-0 | 256X192 | 2 | UNS WAP |
| AX GRE | 2D | GRE | NPW,VBW,EDR, SCIC | 15 | 500 | 20 | 10.42 | 16 | 4-0.5 | 256X192 | 2 | UNS WAP |

Comments: _____

Written: 05/2007
 Revised: 01/2009; 06/2009